



## SP Connect

For use with:

- **SADE SP** tablet/capsule weight sorters
- **SP-4C** tablet/capsule sample weight checker

### Data transfer for improved reporting and accurate analysis

SP Connect is an optional software program that enables you to transfer reports and data from your CI SP unit onto a PC.

With SP Connect enhanced colour reports are automatically saved as secure pdf files for ease of long term storage.

You can also export weight data, including individual weights, into an Excel spreadsheet for further analysis.

### KEY FEATURES

- Reports saved as secure pdf files
- Export data into Excel format
- All files identified with Batch ID
- Enhanced colour graphics

Sample no.	Channel no.	Sample Quantity	Reading no.	Accepted	Reading date & time	Reading value	Reading QC
3	1	0	1	No	25/01/2012 15:27:56	621.10	mg
4	1	0	2	No	25/01/2012 15:28:01	617.60	mg
5	1	0	3	No	25/01/2012 15:28:02	623.10	mg
6	1	0	4	No	25/01/2012 15:28:05	629.40	mg
7	1	0	5	No	25/01/2012 15:28:06	605.40	mg
8	1	0	6	No	25/01/2012 15:28:07	631.30	mg
9	1	0	7	No	25/01/2012 15:28:08	625.00	mg
10	1	0	8	No	25/01/2012 15:28:09	627.40	mg
11	1	0	9	No	25/01/2012 15:28:10	624.00	mg
12	1	0	10	No	25/01/2012 15:28:11	626.00	mg
13	1	0	11	No	25/01/2012 15:28:11	636.10	mg
14	1	0	12	No	25/01/2012 15:28:12	632.20	mg
15	1	0	13	No	25/01/2012 15:28:13	619.30	mg
16	1	0	14	No	25/01/2012 15:28:14	634.30	mg
17	1	0	15	No	25/01/2012 15:28:15	636.40	mg
18	1	0	16	No	25/01/2012 15:28:16	619.60	mg
19	1	0	17	No	25/01/2012 15:28:17	624.60	mg
20	1	0	18	No	25/01/2012 15:28:18	624.20	mg
21	1	0	19	No	25/01/2012 15:28:18	627.70	mg
22	1	0	20	No	25/01/2012 15:28:20	621.50	mg
23	1	0	21	No	25/01/2012 15:28:21	636.10	mg
24	1	0	22	No	25/01/2012 15:28:22	643.00	mg
25	1	0	23	No	25/01/2012 15:28:23	619.60	mg
26	1	0	24	No	25/01/2012 15:28:23	619.50	mg
27	1	0	25	No	25/01/2012 15:28:24	633.40	mg
28	1	0	26	No	25/01/2012 15:28:25	636.50	mg

Reading date & time	Reading
25/01/2012 15:27:56	621.10
25/01/2012 15:28:01	617.60
25/01/2012 15:28:02	623.10
25/01/2012 15:28:05	629.40
25/01/2012 15:28:06	605.40
25/01/2012 15:28:07	631.30
25/01/2012 15:28:08	625.00
25/01/2012 15:28:09	627.40
25/01/2012 15:28:10	624.00
25/01/2012 15:28:11	626.00
25/01/2012 15:28:12	636.10
25/01/2012 15:28:13	632.20
25/01/2012 15:28:14	619.30
25/01/2012 15:28:15	634.30
25/01/2012 15:28:16	636.40
25/01/2012 15:28:17	619.60
25/01/2012 15:28:18	624.60
25/01/2012 15:28:18	624.20
25/01/2012 15:28:18	627.70
25/01/2012 15:28:20	621.50
25/01/2012 15:28:21	636.10
25/01/2012 15:28:22	643.00
25/01/2012 15:28:23	619.60
25/01/2012 15:28:23	619.50
25/01/2012 15:28:24	633.40
25/01/2012 15:28:25	636.50



# SP Connect

Data transfer to PC

## KEY OPERATIONAL BENEFITS

SP Connect enhances the report capabilities of CI's checkweighers.

### Statistical reports automatically saved for archiving

- Sort, Sample and Batch Reports can be transferred
- Calibration reports and system settings can also be saved.

### Security

- Reports are saved as encrypted read only pdf files
- Pdf files protected by 24 character randomly generated password and 40 bit encryption

### Enhanced colour reports

- Graphics presented in colour for ease of interpretation
- Reports can be configured with company name and logo

### Further analysis with Excel

- All report data, including individual weights, can be transferred into a spreadsheet for further analysis

### Connectivity

- SP Connect can collect data from:
  - SADE SP weight sorters
  - SP-4C sample weight checkers
- Multiple CI units can be connected to SP Connect



### Licencing

SP Connect is supplied as a single server licence allowing installation on one PC or network server.

## TECHNICAL DETAILS

Direct Connection	The PC operating SP Connect can be connected directly to one CI unit, using a crossover cable and RJ45 connection.
Router	SP Connect can be set up to run with a router. One or multiple CI units can be connected
Network Server	SP Connect can be set up to run on a network, allowing multiple CI units to be connected via LAN connection.



SP-4C



SADE SP-B40



SADE SP140



SADE SP240



SADE SP440



Brunel Road, Churchfields  
Salisbury, Wiltshire  
SP2 7PX United Kingdom

T +44 (0)1722 424100  
F +44 (0)1722 323222  
E sales@ciprecision.com

[www.ciprecision.com](http://www.ciprecision.com)